

Anti-Muc1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P15941

Item Number: CM0013832



Introduction

Rabbit polyclonal antibody against human MUC1 (Mucin-1, CA15-3) for Western blot, immunohistochemistry (paraffin), and ELISA; cross-reacts with mouse and rat. Ideal for cancer research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	MUC1 Rabbit pAb
Remarks/Alias	EMA, MCD, PEM, PUM, KL-6, MAM6, MCKD, PEMT, CD227, H23AG, MCKD1, MUC-1, ADMCKD, ADTKD2, Ca15-3, ADMCKD1, CA 15-3, MUC-1/X, MUC1/ZD, MUC-1/SEC, MUC1
Species	Human
Gene ID (Human)	4582
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1182-1255 of human MUC1 (P15941)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15941
Protein Weight	122kDa
Shipping	Ice bag