

Anti-Muc16 (Ca125) Polyclonal Antibody For Wb And Elisa - Q8Wxi7

Item Number: CM0020317



Introduction

Rabbit polyclonal antibody targeting human CA125 (MUC16), validated for WB and ELISA. Suitable for ovarian cancer research; detects ~1519 kDa glycoprotein. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	CA125 Rabbit pAb
Remarks/Alias	MUC16; CA125; mucin-16
Species	Human
GeneID (Human)	94025
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 12660-12923 of human CA125 (NP_078966.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8WXI7
Protein Weight	1519kDa
Shipping	Ice bag