

Anti-Muc16/Ca125 Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q8Wxi7

Item Number: CM0006272



Introduction

This rabbit monoclonal antibody against human MUC16/CA125 is validated for WB, IHC-P, IP, and ELISA. Recombinant immunogen ensures specificity, suitable for ovarian cancer research.

[Learn More](#)

Parameter	Value
Product Name	MUC16/CA125 Rabbit mAb
Remarks/Alias	CA125; MUC16
Species	Human
GeneID (Human)	94025
GeneID	94025
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 12660-12923 of human MUC16/CA125 (NP_078966.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human
SWISS	Q8WXI7
Protein Weight	1519kDa
Shipping	Ice bag