

## Anti-Muc2 Monoclonal Antibody For If-P, Ihc-P, Elisa - Q02817

Item Number: CM0004532



### Introduction

Rabbit monoclonal antibody against human MUC2 (Mucin-2), validated for IF-P, IHC-P, and ELISA. Reacts with Human, Mouse, Rat. Immunogen: recombinant fusion protein covering amino acids 858-1259. Suitable for studying intestinal mucosal barrier function.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | MUC2 Rabbit mAb                                                                                                    |
| Remarks/Alias    | MLP; SMUC; MUC-2                                                                                                   |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 4583                                                                                                               |
| GeneID           | 4583                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 858-1259 of human MUC2 (NP_002448.4) |
| Source           | Rabbit                                                                                                             |
| Category         | Monoclonal Antibodies                                                                                              |
| Application      | IF-P, IHC-P, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q02817                                                                                                             |
| Protein Weight   | 551kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                            |