

Anti-Muc2 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q02817

Item Number: CM0012271



Introduction

Rabbit monoclonal antibody against human MUC2, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for intestinal mucin-2 studies in mucosal barrier and immunology research.

[Learn More](#)

Parameter	Value
Product Name	MUC2 Rabbit mAb
Remarks/Alias	MLP, SMUC, MUC-2, MUC2
Species	Human
Gene ID	4583
Immunogen	Synthetic peptide (aa 5000-5100 of human MUC2)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02817
Protein Weight	551 kDa
Shipping	Ice bag