

Anti-Muc2 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q02817

Item Number: CM0015320



Introduction

Rabbit polyclonal antibody targeting human MUC2 (intestinal mucin-2), validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. For mucosal barrier, colon cancer, and copper metabolism research.

[Learn More](#)

Parameter	Value
Product Name	MUC2 Rabbit pAb
Remarks/Alias	MLP; SMUC; MUC-2; MUC2
Species	Human
GeneID (human)	4583
GeneID	4583
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 850-1270 of human MUC2 (NP_002448.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02817
Protein Weight	551 kDa
Shipping	Ice bag