

Anti-Muc4 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q99102

Item Number: CM0010584

Introduction

Rabbit monoclonal antibody against human MUC4, a high-molecular-weight membrane-bound mucin. Validated for WB, IHC-P, and ELISA; cross-reacts with rat. Ideal for oncology and epithelial research.

[Learn More](#)

Parameter	Value
Product Name	MUC4 Rabbit mAb
Remarks/Alias	ASGP; MUC-4; HSA276359; MUC4
Species	Human
GeneID (Human)	4585
GeneID	4585
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 4807-4916 of human MUC4 (Q99102)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q99102
Protein Weight	542kDa
Shipping	Ice bag