

Anti-Muc4 Rabbit Polyclonal Antibody For Wb, Elisa - Q99102

Item Number: CM0024347

Introduction

Rabbit polyclonal antibody against human MUC4, validated for Western blot (WB) and ELISA applications. Targets the membrane-bound mucin involved in epithelial protection and cancer progression. Ideal for oncology research.

[Learn More](#)

Parameter	Value
Product Name	MUC4 Rabbit pAb
Remarks/Alias	ASGP; MUC-4; HSA276359; MUC4
Species	Human
GeneID (Human)	4585
GeneID	4585
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 2100-2169 of human MUC4 (NP_060876.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99102
Protein Weight	542kDa
Shipping	Ice bag