

Anti-Muc6 Polyclonal Antibody For Wb, Elisa - Q6W4X9

Item Number: CM0022022

Introduction

Rabbit polyclonal antibody against human MUC6 (gastric mucin-6). Validated for WB and ELISA; cross-reacts with mouse. Useful for research on mucosal cytoprotection, H. pylori pathology, and tumor marker analysis.

[Learn More](#)

Parameter	Value
Product Name	MUC6 Rabbit pAb
Remarks/Alias	MUC-6; MUC6
Species	Human
Gene ID (Human)	4588
Gene ID	4588
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-298 of human MUC6. (NP_005952.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6W4X9
Protein Weight	257KDa
Shipping	Ice bag