


Anti-Mucin 5Ac Rabbit Polyclonal Antibody For Wb And Elisa - P98088

Item Number: CM0021343



Introduction

Rabbit polyclonal antibody detecting Mucin 5AC, a 586kDa secreted glycoprotein. Validated in WB and ELISA, cross-reacts with human, mouse, rat. Ideal for mucin research, respiratory/gastric studies, and IVD development.

[Learn More](#)

Parameter	Value
Product Name	Mucin 5AC Rabbit pAb
Remarks/Alias	TBM; leB; MUC5; mucin; Mucin 5AC (MUC5AC)
Species	Human
GeneID (Human)	4586
GeneID	4586
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1550-1900 of human Mucin 5AC (MUC5AC) (NP_001291288.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98088
Protein Weight	586kDa
Shipping	Ice bag