

Anti-Multimerin-1 (Mmrn1) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13201

Item Number: CM0029038

Introduction

Rabbit polyclonal antibody targeting human Multimerin-1 (MMRN1), validated for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse MMRN1. Ideal for platelet biology, thrombosis, and vascular hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	MMRN1 Rabbit pAb
Remarks/Alias	ECM; MMRN; GPIa*; EMILIN4; MMRN1
Species	Human
Gene ID	22915
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1083-1228 of human MMRN1 (NP_031377.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13201
Protein Weight	138kDa
Shipping	Ice bag