

Anti-Mum1 Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - Q15306

Item Number: CM0006764



Introduction

Rabbit monoclonal antibody targeting human IRF4 (MUM1/LSIRF), a transcription factor critical for lymphoid cell signaling and immune responses. Validated for IF-P, IHC-P, and ELISA applications.

[Learn More](#)

Parameter	Value
Product Name	MUM1 Rabbit mAb
Remarks/Alias	MUM1, LSIRF, SHEP8, NF-EM5, IRF4
Species	Human
GeneID (Human)	3662
GeneID	3662
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 351-451 of human IRF4 (NP_001182215.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q15306
Protein Weight	52kDa
Shipping	Ice bag