

Anti-Munc 13-1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Upw8

Item Number: CM0023499

Introduction

Rabbit polyclonal antibody targeting human Munc13-1 (UNC13A) for WB and ELISA. Validated in human, mouse, and rat. Ideal for synaptic exocytosis and insulin secretion studies.

[Learn More](#)

Parameter	Value
Product Name	Munc 13-1 Rabbit pAb
Remarks/Alias	Munc13-1; Munc 13-1
Species	Human
GeneID (human)	23025
GeneID	23025
Immunogen	Recombinant protein of human Munc 13-1.
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UPW8
Protein Weight	193kDa
Shipping	Ice bag