

Anti-Mus81 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Ny9

Item Number: CM0027067



Introduction

Rabbit polyclonal antibody targeting human MUS81, a structure-specific endonuclease involved in DNA repair and genome stability. Suitable for WB, IHC-P, ELISA; cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	MUS81 Rabbit pAb
Remarks/Alias	SLX3; MUS81
Species	Human
GeneID (Human)	80198
GeneID	80198
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human MUS81 (NP_079404.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96NY9
Protein Weight	61kDa
Shipping	Ice bag