

Anti-Musashi-1 (Msi1) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O43347

Item Number: CM0004984

Introduction

Musashi-1 (MSI1) rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Reacts with human and mouse. Targets the N-terminus (aa 1-100) of this 39 kDa RNA-binding protein involved in stem cell maintenance and Notch signaling.

[Learn More](#)

Parameter	Value
Product Name	Musashi-1 (MSI1) Rabbit mAb
Remarks/Alias	MSI1; musashi RNA binding protein 1; Musashi-1 (MSI1)
Species	Human
GeneID (human)	4440
GeneID	4440
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MSI1 (O43347).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43347
Protein Weight	39kDa
Shipping	Ice bag