

Anti-Musashi-1 (Msi1) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43347

Item Number: CM0026951

Introduction

Rabbit polyclonal antibody against human Musashi-1 (MSI1), validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for stem cell and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	Musashi-1 (MSI1) Rabbit pAb
Remarks/Alias	MSI1; Musashi-1 (MSI1)
Species	Human
Gene ID	4440
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-362 of human Musashi-1 (MSI1) (NP_002433.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43347
Protein Weight	39kDa
Shipping	Ice bag