


Anti-Musashi-2 (Msi2) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q96Dh6

Item Number: CM0002088



Introduction

High-quality Musashi-2 (MSI2) rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for studying stem cell maintenance and translational regulation.

[Learn More](#)

Parameter	Specification
Product Name	Musashi-2 (MSI2) Rabbit mAb
Remarks/Alias	MSI2H; Musashi-2 (MSI2)
Species	Human
GeneID (human)	124540
GeneID	124540
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Musashi-2 (MSI2) (Q96DH6).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96DH6
Protein Weight	35kDa
Shipping	Ice bag