


Anti-Musashi-2 (Msi2) Rabbit Polyclonal Antibody For Wb, Elisa - Q96Dh6

Item Number: CM0030292



Introduction

High-quality rabbit polyclonal antibody against Musashi-2 (MSI2) for WB and ELISA. Validated for human and mouse, 35kDa protein weight. Ideal for stem cell and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Musashi-2 (MSI2) Rabbit pAb
Remarks/Alias	MSI2H; Musashi-2 (MSI2)
Species	Human
GeneID (Human)	124540
GeneID	124540
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human Musashi-2 (MSI2) (NP_620412.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96DH6
Protein Weight	35kDa
Shipping	Ice bag