


Anti-Muscarinic Achr M2 (Acm2) Monoclonal Antibody For Wb, Ihc-P, Elisa - P08172

Item Number: CM0004021



Introduction

Recombinant rabbit monoclonal antibody against human muscarinic acetylcholine receptor M2 (ACM2). Suitable for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for GPCR signaling and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Muscarinic AChR M2 (ACM2) Rabbit mAb
Remarks/Alias	HM2; Muscarinic AChR M2 (ACM2)
Species	Human
GeneID	1129
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human Muscarinic AChR M2 (ACM2) (P08172)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08172
Protein Weight	52kDa
Shipping	Ice bag