

Anti-Musk Polyclonal Antibody For Wb, Elisa - O15146

Item Number: CM0026360



Introduction

Rabbit polyclonal antibody against human MUSK, validated for WB and ELISA. Cross-reacts with mouse. Ideal for neuromuscular junction research and synaptic signaling studies. 97kDa protein.

[Learn More](#)

Parameter	Value
Product Name	MUSK Rabbit pAb
Remarks/Alias	CMS9; FADS; FADS1; MUSK
Species	Human
GeneID (Human)	4593
GeneID	4593
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 424-783 of human MUSK (NP_001159752.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O15146
Protein Weight	97kDa
Shipping	Ice bag