


Anti-Mut Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P22033

Item Number: CM0028219



Introduction

Rabbit polyclonal antibody against human methylmalonyl-CoA mutase (MUT) for Western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Ideal for mitochondrial metabolism and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	MUT Rabbit pAb
Remarks/Alias	MCM; MUT
Species	Human
GeneID (Human)	4594
Immunogen	Recombinant fusion protein containing amino acids 451-750 of human MUT (NP_000246.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22033
Protein Weight	83kDa
Shipping	Ice bag