

Anti-Mutyh Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q9Uif7

Item Number: CM0023266



Introduction

Rabbit polyclonal antibody against human MUTYH (MYH), a DNA repair enzyme. Validated for IHC-P, IF/ICC, ELISA, and cross-reacts with rat. Ideal for oxidative damage and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MUTYH Rabbit pAb
Remarks/Alias	MYH; MUTYH
Species	Human
Gene ID (Human)	4595
Gene ID	4595
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human MUTYH (NP_001041636.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UIF7
Protein Weight	60kDa
Shipping	Ice bag