


Anti-Mvb12A Polyclonal Antibody For Wb, Elisa - Q96Ey5

Item Number: CM0018986



Introduction

MVB12A Rabbit pAb, a polyclonal antibody raised in rabbit against human MVB12A, is suitable for WB and ELISA. It cross-reacts with rat and is ideal for ESCRT-I complex and endocytic trafficking research.

[Learn More](#)

Parameter	Value
Product Name	MVB12A Rabbit pAb
Remarks/Alias	CFBP; FAM125A; MVB12A
Species	Human
GeneID (human)	93343
GeneID	93343
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-273 of human MVB12A (NP_612410.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q96EY5
Protein Weight	29kDa
Shipping	Ice bag