

Anti-Mvk Polyclonal Antibody For Wb, Elisa - Q03426

Item Number: CM0027403



Introduction

Rabbit polyclonal antibody against human mevalonate kinase (MVK), validated for Western blot and ELISA. Ideal for cholesterol metabolism studies and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	MVK Rabbit pAb
Remarks/Alias	MK; LRBP; MVLK; POROK3; MVK
Species	Human
GeneID (Human)	4598
GeneID	4598
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-300 of human MVK (NP_000422.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q03426
Protein Weight	42kDa
Shipping	Ice bag