


Anti-Mxd1 Rabbit Polyclonal Antibody (Wb, Elisa) - Q05195

Item Number: CM0030600



Introduction

Rabbit polyclonal antibody targeting human MXD1 (MAD). Validated in WB and ELISA applications. Suitable for research on transcriptional repression, MYC antagonism, and telomerase regulation. Cross-reacts with human and rat.

[Learn More](#)

Parameter	Value
Product Name	MXD1 Rabbit pAb
Remarks/Alias	MAD; MAD1; BHLHC58; MXD1
Species	Human
Gene ID (Human)	4084
Gene ID	4084
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human MXD1 (NP_002348.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q05195
Protein Weight	25kDa
Shipping	Ice bag