

Anti-Mxi1 Rabbit Polyclonal Antibody For Wb And Elisa - P50539

Item Number: CM0015834

Introduction

Anti-MXI1 rabbit polyclonal antibody validated for Western blot and ELISA. Detects human, mouse, and rat MXI1, a transcriptional repressor antagonizing MYC. Suitable for cancer and transcription regulation studies.

[Learn More](#)

Parameter	Description
Product Name	MXI1 Rabbit pAb
Remarks/Alias	MXI; MAD2; MXD2; bHLHc11; MXI1
Species	Human
Gene ID (Human)	4601
Gene ID	4601
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-295 of human MXI1 (NP_569157.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50539
Protein Weight	26kDa
Shipping	Ice bag