

Anti-Mybbp1A Rabbit Polyclonal Antibody For Wb - Q9Bqg0

Item Number: CM0026381

Introduction

Rabbit polyclonal antibody targeting MYBBP1A (P160/PAP2). Validated for WB, IHC-P, IF/ICC, ELISA; reacts with human, mouse, rat. Ideal for transcription and ribosome biogenesis research.

[Learn More](#)

Parameter	Value
Product Name	MYBBP1A Rabbit pAb
Remarks/Alias	P160; PAP2; Pol5; MYBBP1A
Species	Human
Gene ID (Human)	10514
Gene ID	10514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1104-1328 of human MYBBP1A (NP_055335.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BQG0
Protein Weight	149kDa
Shipping	Ice bag