

Anti-Mybpc2 Polyclonal Antibody For Wb, Elisa - Q14324

Item Number: CM0017318



Introduction

MYBPC2 rabbit polyclonal antibody for Western blot and ELISA applications. Detects human, mouse, and rat MYBPC2. Useful for investigating muscle protein function and contraction.

[Learn More](#)

Parameter	Value
Product Name	MYBPC2 Rabbit pAb
Remarks/Alias	MYBPC; MYBPCF; fsMyBP-C; MYBPC2
Species	Human
GeneID (human)	4606
GeneID	4606
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 912-1141 of human MYBPC2 (NP_004524.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14324
Protein Weight	128kDa
Shipping	Ice bag