


## Anti-Mybpc2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14324

Item Number: CM0025994



### Introduction

MYBPC2 Rabbit pAb, a polyclonal antibody for WB, IHC-P, ELISA, recognizes myosin-binding protein C fast-type (128kDa) in human, mouse, and rat. Ideal for muscle contraction and structural protein studies.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | MYBPC2 Rabbit pAb                                                                                                    |
| Remarks/Alias    | MYBPC; MYBPCF; fsMyBP-C; MYBPC2                                                                                      |
| Species          | Human                                                                                                                |
| GeneID           | 4606                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 912-1141 of human MYBPC2 (NP_004524.3) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, IHC-P, ELISA                                                                                                     |
| Cross Reactivity | Mouse, Rat                                                                                                           |
| SWISS            | Q14324                                                                                                               |
| Protein Weight   | 128kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                              |