


Anti-Mybpc3 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14896

Item Number: CM0016963



Introduction

Rabbit polyclonal antibody against human MYBPC3 (cardiac myosin-binding protein C), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for cardiac muscle and cardiomyopathy research.

[Learn More](#)

Parameter	Value
Product Name	MYBPC3 Rabbit pAb
Remarks/Alias	FHC; CMH4; CMD1MM; LVNC10; MYBP-C; cMyBP-C; MYBPC3
Species	Human
Gene ID (Human)	4607
Gene ID	4607
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-400 of human MYBPC3 (NP_000247.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14896
Protein Weight	141kDa
Shipping	Ice bag