


Anti-Mycbp Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99417

Item Number: CM0026175



Introduction

Anti-MYCBP rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Targets human c-Myc-binding protein (AMY-1), a key regulator of MYC transcriptional activity. Ideal for cancer pathway studies.

[Learn More](#)

Parameter	Value
Product Name	MYCBP Rabbit pAb
Remarks/Alias	AMY-1; MYCBP
Species	Human
Gene ID (Human)	26292
Gene ID	26292
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-103 of human MYCBP (NP_036465.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q99417
Protein Weight	12kDa
Shipping	Ice bag