


Anti-Myd88 Polyclonal Antibody For Wb, If/Icc, Elisa - Q99836

Item Number: CM0014634



Introduction

Rabbit polyclonal antibody targeting human MyD88 (SWISS: Q99836). Validated for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for TLR/IL-1R innate immune signaling and inflammatory response studies.

[Learn More](#)

Parameter	Value
Product Name	MyD88 Rabbit pAb
Remarks/Alias	WM1, IMD68, MYD88D, MyD88
Species	Human
GeneID (human)	4615
GeneID	4615
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 20-140 of human MyD88 (NP_002459.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99836
Protein Weight	33 kDa
Shipping	Ice bag