


Anti-Myd88 Rabbit Monoclonal Antibody [Ko Validated] For Wb, Ihc-P, Elisa - Q99836

Item Number: CM0003029



Introduction

KO validated rabbit monoclonal antibody against human MyD88, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for innate immunity and TLR signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MyD88 Rabbit mAb
Remarks/Alias	WM1; IMD68; MYD88D; 88
Species	Human
Gene ID (Human)	4615
Gene ID	4615
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human MyD88 (NP_002459.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99836
Protein Weight	15kDa/20kDa/28kDa/31kDa/33kDa/34kDa
Shipping	Ice bag