

Anti-Myd88 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99836

Item Number: CM0015939



Introduction

Rabbit polyclonal antibody to MyD88, validated for WB, IHC-P, ELISA in human, mouse, rat. MyD88 is a key adapter in TLR/IL-1R innate immune signaling pathways. Ideal for immunology and inflammatory disease research.

[Learn More](#)

Parameter	Value
Product Name	MyD88 Rabbit pAb
Remarks/Alias	WM1; IMD68; MYD88D; MyD88
Species	Human
GeneID (human)	4615
GeneID	4615
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-309 of human MyD88 (NP_002459.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99836
Protein Weight	33kDa
Shipping	Ice bag