

Anti-Myd88 Rabbit Polymab® Monoclonal Antibody For Wb, Ip, Elisa - Q99836

Item Number: CM0004886



Introduction

High-quality MyD88 Rabbit PolymAb® monoclonal antibody for WB, IP, ELISA. Detects human and mouse MyD88, a critical TLR/IL-1R innate immune signaling adapter. 33kDa. Q99836.

[Learn More](#)

Parameter	Value
Product Name	MyD88 Rabbit PolymAb®
Remarks/Alias	WM1; IMD68; MYD88D
Species	Human
GeneID	4615
Immunogen	Synthetic peptide corresponding to amino acids 50-150 of human MyD88 (NP_002459.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99836
Protein Weight	33kDa
Shipping	Ice bag