

Anti-Myef2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9P2K5

Item Number: CM0020821



Introduction

High-quality MYEF2 rabbit polyclonal antibody for WB, IHC-P, IF/ICC, and ELISA. Detects MYEF2 in human, mouse, and rat samples. Targets the myelin gene transcriptional repressor, supporting neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	MYEF2 Rabbit pAb
Remarks/Alias	MEF-2, MST156, myEF-2, MSTP156, HsT18564, MYEF2
Species	Human
GeneID	50804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human MYEF2 (NP_001288139.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2K5
Protein Weight	64kDa
Shipping	Ice bag