

Anti-Myelin Basic Protein (Mbp) Monoclonal Antibody For Wb - P02686

Item Number: CM0010485



Introduction

Recombinant rabbit monoclonal antibody targeting human Myelin Basic Protein (MBP). Validated for WB, IF-P, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for neuroscience and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Myelin Basic Protein Rabbit mAb
Remarks/Alias	MBP
Species	Human
GeneID (Human)	4155
GeneID	4155
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02686
Protein Weight	33kDa
Shipping	Ice bag