

Anti-Myelin Basic Protein (Mbp) Rabbit Polyclonal Antibody For Wb, Ihc-P, If-P, Elisa - P02686

Item Number: CM0019669



Introduction

High-quality rabbit polyclonal antibody targeting human Myelin Basic Protein (MBP). Validated for WB, IHC-P, IF-P, and ELISA, with cross-reactivity to mouse and rat. Ideal for neuroscience and demyelination research.

[Learn More](#)

Parameter	Value
Product Name	Myelin Basic Protein Rabbit pAb
Remarks/Alias	MBP; Myelin Basic Protein
Species	Human
GeneID (Human)	4155
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-165 of human Myelin Basic Protein (NP_001020272.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02686
Protein Weight	33 kDa
Shipping	Ice bag