

Anti-Myelin Basic Protein Mouse Monoclonal Antibody For Wb - P02686

Item Number: CM0012780



Introduction

Anti-Myelin Basic Protein (MBP) mouse monoclonal antibody for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Detects a 33 kDa band. Ideal for neurobiology and multiple sclerosis research.

[Learn More](#)

Parameter	Value
Product Name	Myelin Basic Protein Mouse mAb
Remarks/Alias	MBP
Species	Human
GeneID (Human)	4155
GeneID	4155
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-304 of human Myelin Basic Protein (NP_001020272.1).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02686
Protein Weight	33kDa
Shipping	Ice bag