

Anti-Myelin Oligodendrocyte Glycoprotein Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q16653

Item Number: CM0012557

Introduction

Rabbit monoclonal antibody targeting human Myelin Oligodendrocyte Glycoprotein (MOG), a CNS-specific myelin protein and rubella virus receptor. Suitable for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	Myelin oligodendrocyte glycoprotein Rabbit mAb
Remarks/Alias	BTN6; BTNL11; MOGIG2; NRCLP7; Myelin oligodendrocyte glycoprotein
Species	Human
GeneID (Human)	4340
GeneID	4340
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16653
Protein Weight	28kDa
Shipping	Ice bag