

Anti-Myelin Oligodendrocyte Glycoprotein Rabbit Polyclonal Antibody For Wb And Elisa - Q16653

Item Number: CM0030651

Introduction

Rabbit polyclonal antibody against human Myelin Oligodendrocyte Glycoprotein (MOG, Q16653) for WB and ELISA. Cross-reacts with mouse and rat; recombinant immunogen (aa 30-154).

[Learn More](#)

Parameter	Value
Product Name	Myelin oligodendrocyte glycoprotein Rabbit pAb
Remarks/Alias	BTN6; BTNL11; MOGIG2; NRCLP7; Myelin oligodendrocyte glycoprotein
Species	Human
GeneID	4340
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-154 of human Myelin oligodendrocyte glycoprotein (NP_996532.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16653
Protein Weight	28kDa
Shipping	Ice bag