

Anti-Myelin Protein Zero (Mpz) Monoclonal Antibody For Wb - P25189

Item Number: CM0007560



Introduction

Rabbit monoclonal antibody targeting human Myelin Protein Zero (MPZ) for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for peripheral nervous system myelination and neuropathy research.

[Learn More](#)

Parameter	Description
Product Name	Myelin Protein Zero (MPZ) Rabbit mAb
Remarks/Alias	P0; CHM; DSS; MPP; CHN2; CMT1; CMT1B; CMT2I; CMT2J; CMT4E; CMTDI3; CMTDID; HMSNIB; Myelin Protein Zero (MPZ)
Species	Human
GeneID (Human)	4359
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 149-248 of human Myelin Protein Zero (MPZ) (NP_000521.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P25189
Protein Weight	28kDa
Shipping	Ice bag