

Anti-Myelin Protein Zero (Mpz) Polyclonal Antibody For Wb, Elisa - P25189

Item Number: CM0021694



Introduction

Rabbit polyclonal antibody against human Myelin Protein Zero (MPZ) for Western blot (WB) and ELISA. Validated in human, mouse, rat. Ideal for peripheral nerve myelination research.

[Learn More](#)

Parameter	Details
Product Name	Myelin Protein Zero (MPZ) Rabbit pAb
Remarks/Alias	P0; CHM; DSS; MPP; CHN2; CMT1; CMT1B; CMT2I; CMT2J; CMT4E; CMTDI3; CMTDID; HMSNIB; Myelin Protein Zero (MPZ)
Species	Human
Gene ID (Human)	4359
Gene ID	4359
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-153 of human Myelin Protein Zero (MPZ) (NP_000521.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25189
Protein Weight	28kDa
Shipping	Ice bag