


Anti-Myeloblastin (Proteinase 3/Pr3) Rabbit Monoclonal Antibody For If-P, Elisa - P24158

Item Number: CM0003868



Introduction

High-quality rabbit monoclonal antibody targeting human Proteinase 3 (PR3/Myeloblastin). Validated for IF-P and ELISA. Recognizes a neutrophil serine protease, a key c-ANCA autoantigen. - P24158

[Learn More](#)

Parameter	Details
Product Name	Proteinase 3/PR3 Rabbit mAb
Remarks/Alias	MBN; MBT; NP4; P29; PR3; ACPA; AGP7; NP-4; PR-3; CANCA; C-ANCA; Proteinase 3/PR3
Species	Human
GeneID (Human)	5657
GeneID	5657
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Proteinase 3/PR3 (P24158).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Human
SWISS	P24158
Protein Weight	28kDa
Shipping	Ice bag