

Anti-Myeloperoxidase (Mpo) Rabbit Mab For Wb, If/Icc, Ihc-P, Elisa - P05164

Item Number: CM0006460



Introduction

Rabbit monoclonal antibody against human myeloperoxidase (MPO), validated for WB, IF/ICC, IHC-P, and ELISA. Detects MPO from human, mouse, and rat samples—ideal for inflammation and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	Myeloperoxidase (MPO) Rabbit mAb
Remarks/Alias	MPO; myeloperoxidase; Myeloperoxidase (MPO)
Species	Human
GeneID (human)	4353
GeneID	4353
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05164
Protein Weight	84 kDa
Shipping	Ice bag