

Anti-Myeloperoxidase (Mpo) Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P05164

Item Number: CM0013890



Introduction

Rabbit polyclonal antibody against human Myeloperoxidase (MPO), validated for WB, IHC-P, IF-P, and ELISA. Detects MPO in human, mouse, and rat samples. Suitable for studies of neutrophil extracellular trap (NET) formation, innate immunity, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	Myeloperoxidase (MPO) Rabbit pAb
Remarks/Alias	MPO; Myeloperoxidase (MPO)
Species	Human
Gene ID (Human)	4353
Gene ID	4353
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 606-745 of human Myeloperoxidase (MPO) (NP_000241.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05164
Protein Weight	84kDa
Shipping	Ice bag