

Anti-Myh1 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P12882

Item Number: CM0025873

Introduction

High-quality rabbit polyclonal antibody against human MYH1 (myosin-1, MyHC-2x). Validated for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for hearing function and muscle research.

[Learn More](#)

Parameter	Value
Product Name	MYH1 Rabbit pAb
Remarks/Alias	MYHa; HEL71; MYHSA1; MyHC-2x; MyHC-2X/D; MYH1
Species	Human
GeneID (Human)	4619
GeneID	4619
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human MYH1 (NP_005954.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12882
Protein Weight	223 kDa
Shipping	Ice bag