

Anti-Myh13 Polyclonal Antibody For Wb, If-P, Elisa - Q9Ukx3

Item Number: CM0017603

Introduction

Rabbit polyclonal antibody against human MYH13 (Myosin-13), validated for Western blot, immunofluorescence, and ELISA; cross-reacts with mouse and rat. Ideal for studying striated muscle fast-twitch contraction in extraocular and laryngeal muscles.

[Learn More](#)

Parameter	Value
Product Name	MYH13 Rabbit pAb
Remarks/Alias	MyHC-eo; MyHC-IIL; MYH13
Species	Human
GeneID (Human)	8735
GeneID	8735
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1240-1320 of human MYH13 (NP_003793.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UKX3
Protein Weight	224kDa
Shipping	Ice bag