


Anti-Myh14 Polyclonal Antibody For Wb, Elisa - Q7Z406

Item Number: CM0029815



Introduction

Rabbit polyclonal antibody against human MYH14 (myosin-14), suitable for WB and ELISA, cross-reacts with human and rat. A key tool for studying non-muscle myosin IIc in cell shape and cytokinesis.

[Learn More](#)

Parameter	Value
Product Name	MYH14 Rabbit pAb
Remarks/Alias	DFNA4, MHC16, MYH17, PNMHH, DFNA4A, myosin, FP17425, NMHC II-C, NMHC-II-C, MYH14
Species	Human
GeneID (Human)	79784
GeneID	79784
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 600-700 of human MYH14 (NP_079005.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q7Z406
Protein Weight	228kDa
Shipping	Ice bag