

Anti-Myh4 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Y623

Item Number: CM0024560

Introduction

Rabbit polyclonal antibody targeting human MYH4 (Myosin-4, MyHC-2b). Validated for WB, IF-P, IHC-P, ELISA; cross-reacts with human, mouse, rat. For skeletal muscle contraction studies.

[Learn More](#)

Parameter	Value
Product Name	MYH4 Rabbit pAb
Remarks/Alias	MYH2B; MyHC-2B; MyHC-IIb; MYH4
Species	Human
Gene ID (Human)	4622
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1500-1600 of human MYH4 (NP_060003.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y623
Protein Weight	223 kDa
Shipping	Ice bag